

BrewerAccess High-Low EXAM TABLE

Designed to take patients to a new height of care.



Brewer
more possibilities

Durability

- 3-year warranty

Safety

- 18" low height
- Easy on / off for patients
- Safety grab bars for patient stability and positioning
- 450 lb. weight capacity*
- Safety, extra low voltage
- Standard on / off switch

Function

- Patented pass-through side drawers
- Largest storage capacity available
- Pass-through work surface on side of table
- Backrest adjusts for easy patient positioning (pneumatic or power back available)
- Additional paper storage
- Return to chair feature (on power back model)
- Easy to clean, seamless upholstery
- Standard foot control

Patient Comfort

- Easy lateral wheelchair transfers
- Safety grab bars
- Multi-density support upholstery
- Adjustable ergonomic stirrups
- Toe board / leg extension / plastic debris tray
- Premium, plush upholstery available

Aesthetics

- 18 standard upholstery color options
- 6 premium, plush upholstery color options
- Custom colors available

Weight Capacity

- 450 lbs.*

Table Dimensions

- Height: 18" - 37"
- Length: 57" extends to 72"
- Upholstery width: 28"
- Leg extension length: 15"

Covered Paper Roll Storage / Paper Roll Holder

- Storage capacity: 23" x 23" x 3.25"
- Storage for up to 6 additional paper rolls
- Accepts 21" x 3" and 18" x 3" paper rolls

Storage

- Front drawer: 19.7" x 19.2" x 4"
- Side drawers: 22.7" x 20.2" x 7"
- Pass-through side drawers accessible from right or left side of table
- Covered paper roll storage

Certifications

- UL listed, CUL listed

Options

- Power back / return to chair
- Premium plush upholstery
- 5° pelvic tilt / 60-Watt drawer heater in front top drawer

- | | |
|-------|--|
| #6000 | pneumatic back |
| #6001 | pneumatic back with pelvic tilt and drawer heater |
| #6500 | power back, return to chair |
| #6501 | power back, return to chair with pelvic tilt and drawer heater |

accessories

- | | |
|---------|---------------------------------------|
| #100407 | front drawer dividers |
| #98464 | stainless steel treatment pan |
| #100220 | urology drain pan system |
| #99504 | deluxe articulating knee crutches |
| #100230 | deluxe articulating knee crutch clamp |

*under normal use

All Brewer Access High-Low Exam Tables come with a 3-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company
13901 Main Street, Menomonee Falls, WI 53051
1.888.Brewer.1 www.brewercompany.com

ACCESS HIGH-LOW EXAM TABLE

Brewer
more possibilities