

DATE Monday, April 6, 2009

ITEM	3LE300H	DESCRIPTION	Four-section beds
------	---------	-------------	-------------------

SAMPLE IMAGE



- Hydraulic dual-sided mattress platform height adjustment from 43 cm to 80 cm.
- Removable galvanized grid platform surface.
- Crank-gear operated backrest and knee-break adjustment.
- Overall dimensions cm 216x93.



It is not authorized to copying, however in part, or using of the present drawing. The Company reserves the right to introduce all technical and/or aesthetical modifications deemed necessary to their products without prior notice.

Industrie Guido Malvestio S.p.a. - 35010 Villanova (Padova) Italy - Via Caltana, 121 - Tel. 049/9299511 - Fax 049/9299500 - www.malvestio.it
e-mail: info@malvestio.it - cap. soc. 2.450.000,00 € i.v. - c.c.i.a.a. padova 119204 - reg. imprese, c.f. e p.iva 00197370281