

DATE Monday, April 6, 2009

ITEM **3LA860Q**

DESCRIPTION **Electric bed, 4-section, 3 of which jointed and 1 fixed centrally. Mattress platform made of ABS.**

SAMPLE IMAGE



Base frame of welded and epoxy-coated steel, with wheels provided with thermoformed ABS covers. Frame provided with proper housings for fitting patient lifting pole and I.V. poles at the 4 corners. Height adjustment with maximum lifting load of 250 kg. Bed movements push-button control panel (Supervisor) with dedicated software, for the exclusive use of specialized personnel, with inhibition of the controls of the secondary push-button handset in use to the patient. Patient's secondary push-button handset can be hooked to safety sides. Push-button panels provided with on/off leds and sound warning for signalling the inhibited functions. Backrest and knee-break with automatic 16,5 cm translation during their movement. Emergency CPR. Leg rest section adjustable manually by means of a gas spring. Bed stripper and push-button panel holder are integrated and extractable. Removable bed ends in plastic material with coloured inserts. Sound alarm system for warning in case of connection to the power net and wheels not blocked. Epoxy powders coating finish, with antibacterial additive.



It is not authorized to copying, however in part, or using of the present drawing. The Company reserves the right to introduce all technical and/or aesthetical modifications deemed necessary to their products without prior notice.

Industrie Guido Malvestio S.p.a. - 35010 Villanova (Padova) Italy - Via Caltana, 121 - Tel. 049/9299511 - Fax 049/9299500 - www.malvestio.it
e-mail: info@malvestio.it - cap. soc. 2.450.000,00 € i.v. - c.c.i.a.a. padova 119204 - reg. imprese, c.f. e p.iva 00197370281