

DATE Monday, April 6, 2009

ITEM **3LA900Q**

DESCRIPTION **Electric bed, 4-section, 3 of which jointed and 1 fixed centrally. Mattress platform made of ABS.**

SAMPLE IMAGE



Base frame of welded and epoxy-coated steel, with wheels provided with thermoformed ABS covers. Frame provided with proper housings for fitting patient lifting pole and I.V. poles at the 4 corners. Head end mounted onto the base frame so as to protect the wall during trendelenburg and reverse-trendelenburg movements. Height adjustment with maximum lifting load of 250 kg. Bed movements push-button control panel (Supervisor) with dedicated software, for the exclusive use of specialized personnel, with inhibition of the controls of the secondary push-button handset in use to the patient. Patient



It is not authorized to copying, however in part, or using of the present drawing. The Company reserves the right to introduce all technical and/or aesthetical modifications deemed necessary to their products without prior notice.

Industrie Guido Malvestio S.p.a. - 35010 Villanova (Padova) Italy - Via Caltana, 121 - Tel. 049/9299511 - Fax 049/9299500 - www.malvestio.it
e-mail: info@malvestio.it - cap. soc. 2.450.000,00 € i.v. - c.c.i.a.a. padova 119204 - reg. imprese, c.f. e p.iva 00197370281