

DATE Monday, April 6, 2009

ITEM	330600	DESCRIPTION	comodino bifronte
------	--------	-------------	-------------------

SAMPLE IMAGE



comodino bifronte con struttura in lamiera di acciaio verniciata a polveri, piano superiore, frontali cassetti e ante in laminato stratificato con maniglie ad incasso. dotato di portabottiglie esterno posto a lato del paziente, cassetto passante con contenitore estraibile, antine apribili a 1850 su entrambi i lati . alla base ruote in gomma con anello paracolpi, e fondo forato per l'aerazione. dimensioni cm 53x43x72,5h.



It is not authorized to copying, however in part, or using of the present drawing. The Company reserves the right to introduce all technical and/or aesthetical modifications deemed necessary to their products without prior notice.

Industrie Guido Malvestio S.p.a. - 35010 Villanova (Padova) Italy - Via Caltana, 121 - Tel. 049/9299511 - Fax 049/9299500 - www.malvestio.it
 e-mail: info@malvestio.it - cap. soc. 2.450.000,00 € i.v. - c.c.i.a.a. padova 119204 - reg. imprese, c.f. e p.iva 00197370281