

DATE Monday, April 6, 2009

ITEM **331450**

DESCRIPTION **comodino monofronte con frigorifero da 35 lt.  
con sede per carrello servitore.**

SAMPLE IMAGE



comodino con struttura in lamiera d'acciaio verniciato a polveri, con due piani di appoggio dotati di bordo perimetrale di contenimento liquidi, frontale ed anta in laminato stratificato con maniglie antitrauma in materiale plastico morbido. dotato di frigorifero da 35 l interno silenzioso a sbrinamento automatico con antina girevole a 1100, di portabottiglie esterno posto a lato del paziente e di cassetto con contenitore estraibile. alla base 2 ruote in gomma e doppio piedino di appoggio antiscivolo. dotato di sede laterale per l'alloggiamento del tavolino servitore art. 331490 (non compreso). dimensioni cm 55x42x93h



It is not authorized to copying, however in part, or using of the present drawing. The Company reserves the right to introduce all technical and/or aesthetical modifications deemed necessary to their products without prior notice.

**Industrie Guido Malvestio S.p.a.** - 35010 Villanova (Padova) Italy - Via Caltana, 121 - Tel. 049/9299511 - Fax 049/9299500 - www.malvestio.it  
e-mail: info@malvestio.it - cap. soc. 2.450.000,00 € i.v. - c.c.i.a.a. padova 119204 - reg. imprese, c.f. e p.iva 00197370281