

DATE Monday, April 6, 2009

ITEM **330840**

DESCRIPTION **comodino bifronte con piano servitore  
regolabile a leggio e bordo di contenimento**

SAMPLE IMAGE



comodino bifronte con struttura in lamiera di acciaio verniciata a polveri, piano superiore con bordo perimetrale di contenimento liquidi, , frontali cassetti e ante in laminato stratificato con maniglie in materiale plastico antitrauma. dotato di portabottiglie esterno posto a lato del paziente, cassetto passante con contenitore estraibile, antine apribili a 1850 su entrambi i lati e di vassoio servitore in abs, con bordo perimetrale di contenimento liquidi, scorrevole, regolabile in altezza e inclinabile a leggio su entrambi i fronti. alla base ruote in gomma con anello paracolpi, e fondo forato per l'aereazione. dimensioni cm 62x48,5x72,5/110h



It is not authorized to copying, however in part, or using of the present drawing. The Company reserves the right to introduce all technical and/or aesthetical modifications deemed necessary to their products without prior notice.

Industrie Guido Malvestio S.p.a. - 35010 Villanova (Padova) Italy - Via Caltana, 121 - Tel. 049/9299511 - Fax 049/9299500 - www.malvestio.it  
e-mail: info@malvestio.it - cap. soc. 2.450.000,00 € i.v. - c.c.i.a.a. padova 119204 - reg. imprese, c.f. e p.iva 00197370281