

DATE Monday, April 6, 2009

ITEM	DESCRIPTION
384300	poltrona prelievi day-clinic

SAMPLE IMAGE



basamento in tubo rettangolare con 4 ruote e due piedini di bloccaggio con comando a pedale. movimento indipendente di regolazione dello schienale e poggiatesta in acciaio verniciato, comandato da pompa a gas. movimento di trendelenburg ottenuto mediante attuatore elettrico a 24 v. c.c., comandato da pulsantiera posta sul poggiatesta. rivestimento in tessuto plastico lavabile ignifugo tipo similpelle morbida. poggiatesta sagomata, regolabile in altezza ed orientabile su giunto sferico. ingombro massimo cm 170 x 75 dispositivo medico, in classe I, conforme alla direttiva comunitaria 93/42/CEE.



It is not authorized to copying, however in part, or using of the present drawing. The Company reserves the right to introduce all technical and/or aesthetical modifications deemed necessary to their products without prior notice.

Industrie Guido Malvestio S.p.a. - 35010 Villanova (Padova) Italy - Via Caltana, 121 - Tel. 049/9299511 - Fax 049/9299500 - www.malvestio.it
e-mail: info@malvestio.it - cap. soc. 2.450.000,00 € i.v. - c.c.i.a.a. padova 119204 - reg. imprese, c.f. e p.iva 00197370281