

DATE Monday, April 6, 2009

ITEM **345900**

DESCRIPTION **range beds with two columns, 4 sections with translating backrest and electric trendelenburg**

SAMPLE IMAGE



- Electrically operated backrest, knee-break and height adjustment (from 41 cm to 81 cm), with double push-button panel handset and controls lock out.
- Leg-rest with innovative safety micro-switch, to avoid unnatural positions of the knee articulation.
- Removable molded ABS mattress platform.
- Translating backrest designed to reduce pressure in the abdominal region.
- Pair of lateral accessory holder bars.
- Pull-out handset holder with a small bow for blanket-holding use.
- Integrated bed lengthener, with excursion of 18 cm.
- Electrically-operated 14° trendelenburg and reverse trendelenburg with cardiac position.
- Max. safe working load: 250 kg.
- Over-all dimensions: 222x93 cm.



It is not authorized to copying, however in part, or using of the present drawing. The Company reserves the right to introduce all technical and/or aesthetical modifications deemed necessary to their products without prior notice.

Industrie Guido Malvestio S.p.a. - 35010 Villanova (Padova) Italy - Via Caltana, 121 - Tel. 049/9299511 - Fax 049/9299500 - www.malvestio.it
 e-mail: info@malvestio.it - cap. soc. 2.450.000,00 € i.v. - c.c.i.a.a. padova 119204 - reg. imprese, c.f. e p.iva 00197370281